

SIMPLY IMPLANT

TIPOLOGY Less traumatic surgeries for your patients, in comparison with purely conical implants. Simply offers great results both in upper and lower jaw, and it's particularly fitted for D1 and D2 bone classes.

SURGERY Easy and fast insertion due to its shape and precisely calibrated drills, Simply "simplifies" the dentist's job, while adapting to all surgery protocols thanks to a wide range of measures available.

PROSHETICS We offer a wide range of abutments, having taken into consideration every possible situation and shape and every need of the dental technicians, keeping the same prosthetic abutment for all the diameters.

1 SWITCHING PLATFORM

- Expanded contact area between implant surface and abutment
- Enhances the healing process of the gums
- Single prosthetic connection for every diameter with int. hex.

2 MACHINED NECK WITHOUT MICROTHREAD

- Eases implant insertion without the need for crestal tapers
- Reduces the surgical stress on the cortical bone

3 PROGRESSIVE THREAD

- Cutting edge apex, 2nd and 3rd coronal compacting
- Reduces the invasiveness for a better contact with the alveolar bone
- Enhances the osseointegrative response (B.I.C.)

4 ANATOMICAL DESIGN, HYBRID CONTOUR STRUCTURE



- Brand new morphology, perfectly suited to obtain primary and secondary stability
- Load distribution on all the implant body

5 SELF TAPING SYSTEM

- Self taping
- Self drilling

6 ATRAUMATIC PENETRATING APEX

- Allows the implant to penetrate under prepared sockets
- Helps to lift the membrane of the sinus
- Avoids the risk of perforation




**ACID ETCHED
SURFACE**

Simply

Hybrid contour internal hexagon implants

Codes and measurements



| H | Ø 3,5 ■ | Ø 4,0 ■ | Ø 4,5 ■ | Ø 5,0 ■ |
|----------|---|--|--|--|
| 8,5 mm | SIMPLY 3,5-8,5 | SIMPLY 4-8,5 | SIMPLY 4,5-8,5 | SIMPLY 5-8,5 |
| 10,0 mm | SIMPLY 3,5-10 | SIMPLY 4-10 | SIMPLY 4,5-10 | SIMPLY 5-10 |
| 11,5 mm | SIMPLY 3,5-11,5 | SIMPLY 4-11,5 | SIMPLY 4,5-11,5 | SIMPLY 5-11,5 |
| 13,0 mm | SIMPLY 3,5-13 | SIMPLY 4-13 | SIMPLY 4,5-13 | SIMPLY 5-13 |
| 14,5 mm | SIMPLY 3,5-14,5 | SIMPLY 4-14,5 | SIMPLY 4,5-14,5 | - |
| Platform | Normal | Normal | Normal | Normal |

Normal



A version of this implant with a transmucosal machined 1,3 mm collar is available on request



A full machined version of this implant is available on request



Kits and instruments pag. 66
Prosthetic components pag. 91

Surgical protocol Simply implants

Atraumatic osteotomy sequence in compact bone (D1/D2)


Simply implant it's easier to insert than other traditional implants. The initial stability is fast to achieve and a good preparation of the implant site is accomplished through a procedure that takes advantage of the gradual drilling. All the drillings has to be conduct under an abundant external irrigation with sterile solution. Besides, it has to be used a discontinuous drilling technique, in order to avoid the heating of the bone.

***Bone taper required in truly compact bone**


Intervention example with 10mm length implant, positioned at crestal level.

Surgical drills are longer than the implants, the colored insert indicates diameter.


Simply implant \varnothing 3,5 h 10 D1/D2 bone



Simply implant \varnothing 4,0 h 10 D1/D2 bone



Simply implant \varnothing 4,5 h 10 D1/D2 bone



Simply implant \varnothing 5,0 h 10 D1/D2 bone

